



USER MANUAL – EN

**IN 27476 Waist belt with infrared light inSPORTline
Marolino**



SevenSport s.r.o. reserves the right to make any changes and improvements to its product without prior notice. Visit our website www.insportline.eu where you will find the latest version of the manual.

CONTENTS

- SAFETY INSTRUCTIONS..... 3
- SPECIFICATIONS..... 3
- MAINTENANCE AND CLEANING 4
- PRODUCT DESCRIPTION 4
- USE 4
 - POWER SUPPLY 4
 - CONTROLLER 5
- ENVIRONMENT PROTECTION 6
- TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS..... 6

SAFETY INSTRUCTIONS

- Read the manual carefully before first use and keep it for future reference. Follow all warnings and recommendations. Failure to do so exposes you to the risk of injury or damage to the product.
- Use the product only in accordance with the manual. Improper use can result in serious injury. Do not make unauthorized modifications to the product.
- Before starting, it is necessary to put on the belt. Before removing the belt, it is necessary to turn it off.
- Do not use continuously for more than 30 min. You must wait at least 2 hours between uses.
- Do not use the product if your hands are damp or wet. Put the belt on clean clothes or clean and dry skin.
- Do not use on the head or eye area.
- Do not wear jewelry or belts during use.
- The controller and power source must be disconnected before storage and transport.
- Do not cover the belt with a pillow or blanket during use. Heat can build up and damage the product or cause injury.
- Keep away from sources of heat and fire.
- The device is not waterproof. Keep away from water sources. Do not submerge in water.
- Do not plug in with wet or damp hands.
- After use, the device must be turned off and disconnected from the power supply.
- Keep cables and belt away from neck and head.
- For indoor use only. Keep away from water and damp rooms.
- Do not look at the lights unless you are wearing safety glasses.
- Do not use under the influence of drugs, alcohol or medication.
- Do not repair or modify the product yourself. Do not use a damaged product.
- The product is not suitable for pregnant women, people with sensitivity to heat, light, inflammation, rash, eczema or other skin diseases.
- Consult a doctor before use. The device is not a substitute for medical or health care.
- If you experience pain or nausea, stop use immediately and contact a doctor.

SPECIFICATIONS



Number of LEDs	105 pcs
Dimension	130x17 cm
Weight	1 kg
Wavelength	Red, NIR
Rated AC input	100 – 240 V ~ 50/60 Hz
Rated DC power	12V 2A
USB input	5V 3A or 5V 4 A
USB output	12V 1,67A
Power bank (not included)	5V 3A or 5V 4 A

MAINTENANCE AND CLEANING

Wipe with a dry soft cloth. Do not immerse in water or other liquids. Cannot be disinfected with liquid disinfectants.

Store in a dry, shady place away from sunlight and high temperatures.

PRODUCT DESCRIPTION



USE

POWER SUPPLY

WALL ADAPTER

1. Plug the controller into the belt.
2. Plug power source to the controller.
3. Connect the power supply to the electrical current.

USB SOURCE

1. Plug the controller into the belt.
2. Plug power source to the controller.
3. Plug into a 5V 3A or 5V 4A USB power source.


Warning: The belt cannot be powered by a 12V USB adapter.

POWER BANK

The belt can be powered using a 5V power bank. (not included)

1. Plug the controller into the belt.
2. Connect the power bank to the controller.
3. Place the power bank in the pocket.

CONTROLLER

	Wavelength: 660 nm, 850 nm, both types at the same time
	Intensity: L0 – L4
	Timer: 5, 10, 15, 20, 25, 30 min
	Switch: On, 10 Hz / off

Intensity	Power
L0	0,9 W
L1	4,9 W
L2	8,5 W
L3	12,3 W
L4	16,0 W

You can power the product using an electrical source AC 100-240 V, 50/60 Hz or via USB.



ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyards.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions.

Don't put the batteries among house waste but hand them in to the recycling place.

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated here under determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Strakonická street 1151/2c, Prague 150 00, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

"The Buyer who is the End Customer" or simply the "End Customer" is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

"The Buyer who is not the End Customer" is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ
SEVEN SPORT s.r.o.

Registered Office: Strakonická 1151/2c, Praha 5, 150 00, ČR
Headquarters: Dělnická 957, Vítkov, 749 01
Warranty & Service: Čermenská 486, Vítkov 749 01

CRN: 26847264
VAT ID: CZ26847264
Phone: +420 556 300 970
E-mail: eshop@insportline.cz
reklamace@insportline.cz
servis@insportline.cz

Web: www.inSPORTline.cz

SK
inSPORTline s.r.o.

Headquarters, warranty & service center: Električná 6471,
Trenčín 911 01, SK

CRN: 36311723
VAT ID: SK2020177082
Phone: +421(0)326 526 701
E-mail: objednavky@insportline.sk
reklamacie@insportline.sk
servis@insportline.sk

Web: www.inSPORTline.sk

About shipping

